

Hygopac Handover Protocol

This protocol confirms the qualified handover and instruction for the Hygopac rotary sealing device. It represents a constituent part of the installation qualification (IQ). It is recommended that you have this carried out by a qualified adviser from a specialist dealer, who will draw your attention to the special features of the system. Please complete the following document and return the original to Dürr Dental.

Hygopac serial number:

- ☐ Unpacking the device with an inspection for damage and complete content
- ☐ Installation and setup of the device
- ☐ Instruction in the operation of the device
- ☐ Setting the sealing temperature in accordance with the installation and operating instructions
- ☐ Functional test, e.g. with Hygopac Sealcheck/Hygofol and visual inspection of the sealing seam
- ☐ Explanation of the daily/annual validation process (operational qualification (OQ)/performance qualification (PQ))
- ☐ Please provide further information as necessary to complete the protocol. Original to Dürr Dental, 1. duplicate to the specialist dealer, 2. duplicate to be kept by operator

Practice stamp

Supplier stamp

I confirm the handover according to the specifications mentioned above.

Operator name in BLOCK LETTERS

Further operators in BLOCK LETTERS

Further operators in BLOCK LETTERS

Name of the expert advisor in BLOCK LETTERS

Date/signature of operator

Date/signature of operator

Date/signature further operators

Date/Signature of the adviser

DÜRR DENTAL SE
Höpfigheimer Str. 17
74321 Bietigheim-Bissingen
Germany
www.duerrdental.com
info@duerrdental.com

